

# Breakthrough Antibody Test for COVID-19

Detects all 3 antibodies (IgA, IgM, IgG) in just 15 minutes via oral swab

**CovAb™** is an oral sample COVID-19 rapid antibody test. Unlike other antibody tests that require a blood draw, CovAb™ is administered with a simple mouth swab, making it easy and painless for patients of all ages. The kit does not require instruments, and the test is completed on-site with results in just 15 minutes.

CovAb™ stores at room temperature and is manufactured in the U.S.A. Available only through select distributors. **To learn more about becoming a distributor, visit [CovAb.com](https://CovAb.com)**



RESULTS



1. Swab
2. Swish
3. Start timer
4. See results



- Rapid and easy 15-minute test with oral swab
- All-inclusive kit; no instrument required
- 97.6% sensitivity / 98.8% specificity\*
- Test kits stored at room temperature
- Made in the U.S.A.

\*≥15 days from onset of symptoms.



**CovAb.com** |  
**1.877.748.9355**



RD-RatioDiagnostics GmbH  
Westerbachstraße 47  
60489 Frankfurt / Main  
Germany

1. This product has not been FDA cleared or approved, but has been authorized for emergency use by FDA under an EUA for use by authorized laboratories;  
2. This product has been authorized only for detecting the presence of total antibodies to SARS-CoV-2, not for any other viruses or pathogens; and  
3. The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.